

Standard Operating Procedure for trace gas measurements by Proton Transfer Reaction-Mass Spectrometry

Authors

S. Dusanter	IMT Nord Europe
R. Holzinger	Utrecht Univ.
F. Klein	Deutscher Wetterdienst
<u>T. Salameh</u>	IMT Nord Europe
M. Jamar	IMT Nord Europe

Contributions from


M. Muller	IONICON Analytik GmbH
C. Amelynck & B. Verreyken	Royal Belgian Institute for Space Aeronomy
A. Munöz	CEAM
I. Ylivinkka	Univ. of Helsinki

SOP / Outline of presentation

These guidelines contain the following topics:


1.	General introduction	2
2.	Principle of the PTR-MS technique	3
3.	Quality Assurance	7
a.	Blank measurements	7
b.	Calibration	8
c.	Figures of merit	11
4.	Field operation	14
a.	Sampling	15
b.	Operating conditions	16
c.	Checking for proper field operation	17
5.	Data extraction: Retrieval of ambient VMR from mass spectra & reporting	17
6.	References	19
	Annex 1: Non-exhaustive list of compounds detected at specific m/z values	21
	Annex 2: Equilibration time required for zeroing the instrument	22
	Annex 3: NPL gas standard	23
	Annex 4: Examples of ion transmission curves	24
	Annex 5: Examples of humidity-dependent sensitivities	25
	Annex 6: Detector voltage optimization	26
	Annex 7: Evaluation of primary ion purity and distribution	27

PTRMS models covered in this SOP



- **Ion Source:** H_3O^+ & traces of NO^+ , O_2^+

- **Drift Tube:**



- **Ion Transfer Unit:** electrostatic lenses, focusing ion technologies
- **Mass Spectrometer:** QMS, ToFMS

This SOP covers only:

- ***PTR-MS equipped with conventional drift tubes*** - simple reaction kinetics and energetics
- ***Proton transfer using H_3O^+*** - use of other reagent ions such as O_2^+ and NO^+ not included

Quality Assurance

Blank measurements

$$C_{RH^+} = C_{RH^+}^{ambient} - C_{RH^+}^{blank}$$

Zero air: VOC free air produced from ambient air using a catalyst @ 350-400°C; a VOC scrubber containing a solid sorbent, or (iii) a cylinder of dry zero air (high purity, 5.0) with a humidification system containing high purity water

Frequency: hourly blanks recommended. At least one blank every six hours to ensure that diurnal variations in instrumental blanks are well captured

Duration: several tens of minutes / check for blank stability

Ref blank values: blanks larger than reported in the table should be investigated

3. Quality Assurance
 - a. Blank measurements
 - b. Calibration
 - c. Figures of merit

Table 3: Blank values for PTR-ToFMS

Species	Formula	m/z	Blank values (pptv)*
Methanol	CH ₃ OH	33.033	30-90
Acetonitrile	C ₂ H ₃ N	42.034	30-110
Acetaldehyde	C ₂ H ₄ O	45.033	40-270
Acetone	C ₃ H ₆ O	59.049	60-500
Isoprene	C ₅ H ₈	69.070	20-70
Methyl Ethyl Ketone	C ₄ H ₈ O	73.065	10-70
Benzene	C ₆ H ₆	79.054	20-80
Toluene	C ₇ H ₈	93.070	20-140
C8-aromatics	C ₈ H ₁₀	107.086	10-110
C9-aromatics	C ₉ H ₁₂	121.101	5-200
Monoterpenes	C ₁₀ H ₁₆	81.070 + 137.132	60-100 + 10-50

*Range of values observed on several PTR-ToFMS (m/Δm>3000). Reported values are for systems which have been running for at least several days. PTR-QMS may exhibit larger values due to the detection of isobaric species generated in the ion source and desorbing from the inlet material.

Quality Assurance

Calibration

Provide procedures & recommendations

3. Quality Assurance
 - a. Blank measurements
 - b. Calibration
 - c. Figures of merit

zero air and calibration gas (few ppb) sampled until a stable plateau is reached for each VOC signal
Frequency: once a week (campaign < 2 months) / twice a month (> 2 months)

Kinetic approach

$$VMR = \frac{C_{RH^+} \times IF_{RH^+} \times FF_{RH^+}}{C_{H_3O^+}} \left/ \left(\frac{T_{RH^+}}{T_{H_3O^+}} \times \frac{d \times k}{K_0 \times N_0} \times \frac{N^2}{E} \right) \right.$$

- **Use standard mixture** of VOCs w/ known proton transfer rate constants covering the whole mass range interest for ambient measurements
- **Calibrate ion transmission** at low humidity (e.g. 30% RH)
- **Evaluate whether humidity-dependent responses** are observed by varying humidity over 30-90% RH
- **Characterize fragmentation patterns** of species of interest $\rightarrow FF_{RH^+}$

Calibration approach

$$VMR = \frac{C_{RH^+} \times 10^6}{C_{H_3O^+}} \left/ S_N(RH^+) \right.$$

- **Use standard mixture** of all VOCs targeted for ambient measurements
- **Calibrate $S_N(RH^+)$** for all VOCs
- **Evaluate whether humidity-dependent responses** are observed and, if so, apply the Following procedure for correction

$$VMR = \frac{\frac{C_{RH^+}}{C_{H_3O^+} + Xr \times C_{H_3O^+(H_2O)}} \times 10^6}{S_N(RH^+)}$$

Quality Assurance Calibration

Calibration standard from NPL
(Worton et al., AMT, 2023)

- SI-traceable
- 20 compounds tailored to the requirements of PTR-MS
- Expanded uncertainty: 5.1-10.8%

Component	Amount fraction / (μmol/mol)	Component	Amount fraction / (μmol/mol)
Acetaldehyde	1.35 ± 0.14	1,2,4-trifluorobenzene	0.928 ± 0.047
Methanol	1.69 ± 0.17	Toluene	0.967 ± 0.049
Ethanol	1.11 ± 0.12	Hexamethylcyclotrisiloxane	0.98 ± 0.10
Isoprene	1.13 ± 0.06	<i>m</i> -xylene	0.963 ± 0.049
Acetone	1.07 ± 0.06	Octamethylcyclotetrasiloxane	1.11 ± 0.12
Dimethyl sulfide	1.05 ± 0.06	1,2,4-trimethylbenzene	1.02 ± 0.11
Acetonitrile	1.29 ± 0.13	+3-carene	0.92 ± 0.10
Perfluorotributylamine	0.92 ± 0.10	Decamethylcyclopentasiloxane	1.11 ± 0.12
3-buten-2-one	1.14 ± 0.12	1,2,4-trichlorobenzene	1.04 ± 0.11
Butan-2-one	1.16 ± 0.06	Nitrogen	Balance
Benzene	0.845 ± 0.043	-	-

Table A2.1. Composition of the NPL standard.

Compound	Formula	Protonated mass(es)	Fragments
Nitrogen (balance)	N ₂	-	-
Methanol	CH ₄ O	33.033	-
Acetonitrile	CH ₃ CN	42.034	-
Acetaldehyde	C ₂ H ₄ O	45.033	-
Ethanol	C ₂ H ₆ O	47.050	29.039
Acetone	C ₃ H ₆ O	59.049	
Dimethyl sulfide	C ₂ H ₆ S	63.030	
Isoprene	C ₅ H ₈	69.070	41.039
Methylvinylketone (3-buten-2-one)	C ₄ H ₆ O	71.049	
Methylethylketone (Butan-2-one)	C ₄ H ₈ O	73.065	
Benzene	C ₆ H ₆	79.054	
Toluene	C ₇ H ₈	93.070	
<i>m</i> -Xylene	C ₈ H ₁₀	107.086	
1,2,4-Trimethylbenzene	C ₉ H ₁₂	121.101	
1,2,4-Trifluorobenzene	C ₆ H ₃ F ₃	133.026	
3-Carene	C ₁₀ H ₁₆	137.132,	81.070
1,2,4-Trichlorobenzene	C ₆ H ₃ Cl ₃	180.937, 182.934, 184.931	
Hexamethylcyclotrisiloxane (D3)	C ₆ H ₁₈ Si ₃ O ₃	223.064, 224.063, 225.061	207.032, 208.032, 209.029
Octamethylcyclotetrasiloxane (D4)	C ₈ H ₂₄ Si ₄ O ₄	297.083, 298.082, 299.079	281.051, 282.051, 283.048
Decamethylcyclopentasiloxane (D5)	C ₁₀ H ₃₀ Si ₅ O ₅	371.101, 372.101, 372.105, 373.098	355.070, 356.070, 357.067
Perfluorotributylamine	C ₁₂ F ₂₇ N	671.968	413.977

Quality Assurance

Figures of merit

3. Quality Assurance
 - a. Blank measurements
 - b. Calibration
 - c. Figures of merit

Provide methodologies to compute figures of merit w/ various examples

Mass resolution & accuracy

$$R_{FWHM} = \frac{(m/z)_{peak}}{\Delta(m/z)_{FWHM}}$$

Limit of detection (3σ)

$$LOD = 3 \times \frac{\sqrt{C_{AH^+}^{blank}}}{S_m(AH^+)}$$

Meas. precision (1σ)

$$\sigma_{precision} = \sqrt{C_{AH^+}^{ambient} + C_{AH^+}^{blank}}$$

Where $C_{AH^+}^{ambient}$ and $C_{AH^+}^{blank}$ are expressed in counts.

$$RSD(precision) = \frac{\sigma_{precision}}{C_{AH^+}^{ambient} - C_{AH^+}^{blank}}$$

Meas. accuracy (1σ)

$$RSD(accuracy) = \sqrt{\left(\frac{\sigma_{CAL}}{VMR_{CAL}}\right)^2 + \frac{1}{(F_{CAL} + F_{AIR})^2} \left(\frac{F_{AIR}^2}{F_{CAL}^2} \times \sigma_{F_{CAL}}^2 + \sigma_{F_{AIR}}^2 \right)}$$

Field Operation

Sampling

Provide recommendations on various aspects

Height: 2-3 different heights should be tested for a significant period of time

Material: coated stainless steel (Silcosteel or SilcoNert®1000 and Sulfinert or SilcoNert®2000) or polymer after a sufficient passivation time before use (PFA, PTFE)

Sampling line heating: 40-50°C

Flow rate: a few L/min to a few tens of L/min (residence time < a few seconds), turbulent flow preferred

Particle filter: mesh size < 5 µm unless low volatility compounds are targeted

4. Field operation
 - a. Sampling
 - b. Operating conditions
 - c. Checking for proper field operation




Figure 1:
PTR-MS setup for
field operation

Field Operation

Operating conditions

Provide settings for operating conditions & requirements for reagent ions

4. Field operation
 - a. Sampling
 - b. Operating conditions
 - c. Checking for proper field operation

Table 4: Field operating conditions

Parameter	Range	Comments
Temperature (°C)	50-60 (up to 80 if IVOC & SVOC are targeted)	This temperature should be kept higher than that of the PTRMS-inlet to prevent condensation. Higher temperature: less memory effects and improved response time, less losses of semi-volatile VOCs.
Pressure (mbar)	1-4	Higher pressure reduces the diffusion of sampled air into the ion source in certain PTR-MS models, resulting in lower impurity levels of O_2^+ and NO^+
E/N (Td)	120-140	A higher E/N reduces the abundance of $H_3O^+(H_2O)$ but increases the fragmentation of protonated VOCs
PTR-MS inlet Sampling flow rate (SCCM)	50-500	A higher flow rate reduces the impact of instrument related contamination on blank signals
Sampling line flow rate	See section 4.a	A higher flow rate reduces the impact of sampling line related contamination on blank signals

Table 5: Criteria for optimum tuning of the ion source/drift tube

Ion source	
Primary ions purity	$[O_2^+ + NO^+] < 3\%$ of $[H_3O^+ + H_3O^+(H_2O)]$
Primary ions distribution	$[H_3O^+(H_2O)] < 20\%$ of $[H_3O^+ + H_3O^+(H_2O)]$

Field Operation

Checking for proper field operation

Build control charts of relevant metrics to ensure that instruments are running well

4. Field operation
 - a. Sampling
 - b. Operating conditions
 - c. Checking for proper field operation

Table 6: Control metrics for proper field measurements


Metric	Purpose	Criteria
H ₃ O ⁺ ion signal at m/z 21.022	Check cleanliness of reagent ion source and/or a change in detection sensitivity	Not lower than 70% of initial value at installation
[O ₂ ⁺ +NO ⁺] / [H ₃ O ⁺ +H ₃ O ⁺ (H ₂ O)]*	Check for reagent ions purity	As defined in Table 5
[H ₃ O ⁺ (H ₂ O)] / [H ₃ O ⁺ +H ₃ O ⁺ (H ₂ O)]*	Check for reagent ions distribution	As defined in Table 5
R _{FWHM} (section 3c) & Peak shape	Check ToFMS tuning	R _{FWHM} not lower than 70% of initial value at installation
Blank values (sections 3a & 4b)	Check for leaks and/or contamination	Lower than 1.3 times initial values at installation
N ₂ H ⁺ ion signal at m/z 29.014	Check collision energetics in transfer region - N ₂ H ⁺ and CO ₂ H ⁺ only formed via endothermic reactions in transfer system	Minimize at selected ion-source/drift operating conditions
CO ₂ H ⁺ ion signal at m/z 44.998		
Fragmentation reference species** Ex: Diiodobenzene - ratio of m/z 203.943 over m/z 330.848	Check for abnormal ion fragmentation	Deviation lower than 30% from initial value at installation

* Ratios of transmission corrected signals. **Reference compounds continuously added to the drift tube such as diiodobezene in the PerMaScale feature of IONICON instruments


Data extraction: retrieval of ambient VMR & reporting

Provide a general ‘step by step’ procedure to process PTR-MS raw data

Software tools provided with PTR-MS (PTR-MS viewer or Ionicon Data Analyzer from IONICON Analytic GmbH, Tofware from TOFWERK) or by the scientific community (Toftool, PTRwid)



Grauss et al.
(2010)



Kinetic approach

$$VMR = \frac{C_{RH^+} \times IF_{RH^+} \times FF_{RH^+}}{C_{H_3O^+}} \bigg/ \left(\frac{T_{RH^+}}{T_{H_3O^+}} \times \frac{d \times k}{K_0 \times N_0} \times \frac{N^2}{E} \right)$$

Calibration approach

$$VMR = \frac{C_{RH^+}}{C_{H_3O^+}} \times 10^6 \bigg/ S_N(RH^+)$$

Table 7: Reference peaks for mass calibration

Species	Formula	m/z (Da)	Comments
Hydronium	H ₃ ¹⁸ O ⁺	21.022	
Hydronium water-cluster	H ₃ ¹⁶ O ⁺ (H ₂ ¹⁶ O)	37.028	May be saturated
Acetone.H ⁺	C ₃ H ₆ O.H ⁺	59.049	May be missing when blank measurements are performed
Diiodobenzene.H ⁺ and fragment	C ₆ H ₄ I ₂ .H ⁺	330.848	PerMaScale feature from IONICON instruments
	fragment	203.943	

Check EBAS website for templates and vocabulary updates

Additional PTRMS related activities

- Physical dedicated hands-on training to organize for small groups, or during upcoming intercomparisons with other CiGas units
- Possibility to organize intercomparisons (last one dedicated to VOC was for formaldehyde) & in 2019 for PTRMS within EUROCHAMP
- Operational for evaluation of working standards and target gas, in order to link the VOC measurement by PTRMS to primary standards and ensure the state-of-the art traceability to international scales
- Round-robin: send out a NPL PTRMS calibration standard, a dilution system, and/or an unknown target gas
- Audit and Labelling (step 1a & 1b)
- SOP for instruments equipped with new generation drift tubes ??